COLPHARMA

Thermo Green

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Made in P.R.C.

Wuxi Medical Instrument
Factory - Jiangsu, China

EC REP

Lotus Global Co., Ltd
1 Four Seasons Terrace,
West Drayton, Middlesex,
London, UB796G, United
Kingdom.

CE 1023

Thermo Green

Ecological gallium thermometer

EN

IMPORTANT PRECAUTIONS

- This thermometer is intended exclusively for measuring body temperature!
- Do not drop the thermometer and do not expose it to knocks.
- Do not bend, damage or bite the thermometer tip.
- Thermometers with damaged glass must not be used (risk of injury).
- The measurement liquid contained in the thermometer, a gallium alloy, must never come into contact with metals due to a reaction hazard.
- When using the thermometer, take all necessary precautions to prevent the risk of injury.
- Keep out of the reach of children.
- Children or disabled individuals must always use this device under adult supervision..
- Using the thermometer does not replace medical advice. A persistent fever, particularly in children, must be medically treated.

Warning: if you compare temperatures between two thermometers, you may have two different readings even if both thermometers are within their accuracy tolerances. Differences occur in each measurement.

Read the instructions carefully before using this device.

Store in a dry place.

USING THE THERMOMETER

1. Before measuring body temperature, make sure that the gallium alloy contained in the measuring column is between the first and second numbered mark (that is between 35 °C and 36 °C).
2. Clean/disinfect the thermometer before and after use. Follow the instructions in the “Cleaning and Disinfection” section.
3. Position the thermometer according to the type of measurement chosen (armpit, rectum or mouth).
4. It takes at least 5 minutes after positioning for the thermometer to complete the body temperature measurement. To facilitate the temperature reading, move the thermometer slowly back and forth. The level of gallium alloy in the measuring column will indicate the body temperature reached on the temperature scale.

Please note: Immediately after each measurement, the level of gallium alloy in the measuring column must be returned to between 35 °C and 36 °C. To do this, use the thermometer cover as described in “How to reset the thermometer” section.

HOW TO RESET THE THERMOMETER

- Place the thermometer in its shockproof plastic container.
- Screw the cover in tightly.
- Grip the flexible tab with your thumb and forefinger.
- Hold the thermometer downward.
- Shake the thermometer from right to left, moving only the wrist, until the measurement reads between 35 °C and 36 °C.

Please note: If the storage temperature is lower than 20 °C, keep the thermometer at a temperature ranging from 20 °C to 30 °C for at least 1-2 hours or keep the thermometer at a temperature ranging from 30 °C to 35 °C for ten minutes before using it.

MEASUREMENT TYPES/NORMAL BODY TEMPERATURE

Oral - about 37 °C

Physicians prefer to measure temperature from the mouth since it provides fast and accurate results. Insert the thermometer in the mouth

left or right side, near to the root of the tongue. The tip of the thermometer must be in contact with the tissue. Close the mouth and breathe only through the nose to avoid affecting the measurement.

Approximate measurement time: 5 minutes.

Rectal - about 37.6 °C

This is the safest method and particularly suitable for infants and children. Carefully insert the tip of the thermometer into the anus up to a maximum of 2 cm. To facilitate insertion, we recommend lubricating the tip of the thermometer with a little ointment or oil for children's massage; do not use Vaseline. Stop immediately if you encounter any resistance.

Approximate measurement time: 5 minutes.

Axillary - about 36.7 °C.

For more reliable results, oral or rectal temperature measurement is recommended. Dry the armpit before measuring. Place the tip of the thermometer under the armpit and hold the arm tight to the torso.

Approximate measurement time: 5 minutes.

CLEANING AND DISINFECTION

The thermometer has an airtight glass structure and can be cleaned/disinfected without limits. For cleaning, before and after each use, wash the thermometer with cold water and soap. For disinfection, use a cotton pad soaked in 70% of Isopropyl alcohol and disinfect the thermometer. Alternatively, immerse the thermometer in a disinfectant solution (follow the disinfectant manufacturer's instructions). **At the hospital, the thermometer must be disinfected by trained personnel.**

Warning: the temperature of the cleaning and disinfectant solution must not exceed 42 °C.

If the thermometer breaks, the spilled liquid can be collected and removed using a brush, cellulose cloth or cloth moistened with soapy water. Everything can be disposed of safely as household waste.

TECHNICAL SPECIFICATIONS

All thermometers were inspected by the manufacturer to comply with the prescribed measurement tolerances (+0.1 / -0.15 °C).

Type:	glass thermometer with measuring column containing gallium alloy.
Measuring range:	from 35.0 °C to 42.0 °C
Measurement tolerance:	0.1 °C
Measurement accuracy:	+0.1 °C / -0.15 °C
Operating conditions:	from 20 °C to 42 °C
Storage conditions:	-10° C - 42° C

This device complies with the 93/42/EEC Medical Products Directive. EN12470-1 clinical thermometers standard, part 1 (glass thermometers with metallic liquid expansion with maximum operating device).

We reserve the right to apply technical changes.

WARRANTY

This device has a 2-year warranty from the date of purchase. The warranty is valid only with the original purchase receipt Opening or tampering with the device voids the warranty. The warranty does not cover damage caused by improper use, accidents or failure to abide by the instructions for use.

Colpharma s.r.l.
Via A. M. Vicenzi 19/A - 43124 Parma - Italy
info@colpharma.com - www.colpharma.com

Made in P.R.C.
Wuxi Medical Instrument
Factory - Jiangsu, China

EC REP Lotus Global Co., Ltd
Four Seasons Terrace,
West Drayton, Middlesex,
London, UB796G, United
Kingdom.

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