COLPHARMA

Air 110

F A M I L Y

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Air 110

F A M I L Y

Compressed-air aerosol therapy system

EN

- | | | | |
|---|------------------|----|---|
| 1 | Compressor | 8 | Mouthpiece |
| 2 | ON/OFF button | 9 | EASY 2 Nebulizer |
| 3 | Air tube outlet | 10 | Soft nosepiece |
| 4 | Nebulizer holder | 11 | Air tube |
| 5 | Power cable | 12 | Five air filters |
| 6 | Adult mask | 13 | Components and assembly of EASY 2 nebulizer |
| 7 | Paediatric mask | 14 | Air filter replacement |

Please read the instructions carefully before using this device.

BF type part fitted.

Dear Customer,

This device is a high-quality aerosol therapy system, ideal for occasional home use.

It is fitted with the EASY 2 nebulizer, which has been designed to treat respiratory illnesses quickly and effectively.

It can be used with the most common aerosol-therapy medicines. If you have any questions or problems, or if you need to order spare parts, please contact your dealer or the COLPHARMA customer care service.

Alternatively, you can visit the www.colpharma.com website where you will find additional information on all our products.

INTENDED USE

This product has been designed for the inhalation of medicines to treat respiratory illnesses. It can be used by adult and paediatric patients.

Intended users

- Certified medical professionals (physicians, nurses and therapists).
- Patients at home with home healthcare support for treatment under the direction of a qualified physician.

The user must understand the general operation of the device and the contents of the instruction manual.

It is forbidden to use the device on patients who are unconscious or who cannot breathe on their own.

Intended patients

People who suffer from respiratory illnesses and who need to take inhalation medicines.
Suitable both for paediatric and adult use.

Environment

This product has been designed for home use, as well as for medical clinic settings and hospital use.

Side effects

None.

Minimum product lifespan

The minimum product lifespan periods are reported below and are based on the assumption that the product is used to nebulize a saline solution (2 ml) twice a day, for 7.5 minutes at a time, at an ambient temperature of 23 °C.

The lifespan period may vary depending on the usage environment.

- Compressor (main unit): 2 years.
- Nebulizer - air tube - mouthpiece - nosepiece - paediatric mask – adult mask: 1 year.
Frequent product use may shorten the expected lifespan.
- Air filter: 60 days.

Warnings for use

Follow the warnings and precautions reported in the instruction manual.

IMPORTANT PRECAUTIONS

- Follow the instructions for use. This manual contains important information on the operation and safety of this device. Please read the instructions carefully before using the device, and keep them handy for future reference.
- This device must be used exclusively as described in this manual. The manufacturer shall not be liable for any damage caused by improper use.
- The device is not suitable for anaesthesia or pulmonary ventilation.
- The device must only be used with its original accessories as listed in the user manual.
- Do not use the device if you note any damage or anomaly.
- Never open the device.
- This device has been manufactured using sensitive components and it must be handled with care. Please follow the storage and operation instructions detailed in the “Technical Specifications” chapter.
- The device must be protected from:
 - water and humidity
 - extreme temperatures
 - impacts and falls
 - contamination and dust
 - direct sunlight
 - heat and cold
- Only use this device with the medicines prescribed by the physician and follow dosage, duration and frequency instructions for the prescribed treatment.
- Do not bend the nebulizer beyond 45°.
- Do not use the device close to strong magnetic fields such as mobile phones or radio installations. Always maintain a minimum distance of 3.3 m from other equipment when using this device.
- Never use the device together with other medical devices or electrical equipment.
- Never use the device at a temperature above 40° C.
- Store the device and its accessories in a clean place.
- Do not use or store this device if smoke or gas is present.
- Do not use any medicines left in the nebulizer. Always use a fresh new dose.
- Always make sure that nothing is clogging the filter.
- Materials in contact with the human body are tested for biological safety (ISO 10993) and meet the requirements for irritation and sensitisation.

Children must always use this device under adult supervision. Some parts are small and may represent a choking hazard if swallowed. Cables and tubes may represent a choking hazard.

Never make changes to the device unless authorised to do so by the manufacturer. Never disassemble or try to repair the device or its parts without authorisation and never perform any maintenance interventions while the device is being used by a patient.

Electronic devices must be disposed of according to local regulations and never as ordinary household waste.

The use of this device does not replace the medical advice of a physician.

HOW TO PREPARE AND USE THE DEVICE

When using the device for the first time, make sure to clean all parts as described in the section “Cleaning and Disinfection”.

1. Open the nebulizer (9) and pour in the amount of medicine as prescribed by the physician. Make sure not to pass the maximum level mark (Fig. 2/A).
2. Choose the most suitable inhalation accessory, as prescribed by the physician (Fig. 2/B).
3. Connect the nebulizer to the compressor using the air tube (Fig. 2/C) and plug the power cable into an electric outlet, making sure that the supply voltage is in compliance with that stated in this manual or on the device nameplate.
4. Press the ON/OFF (2) button to turn on the device. Place the mouthpiece in the mouth (8), let the nosepiece rest (10) under the nostrils or the mask (6 or 7) on the face, making sure to cover mouth and nose.
5. Sit in a relaxed position with the upper body straight. Do not lie down while inhaling.
6. After completing the treatment session as recommended by the physician, press the ON/OFF (2) button to turn off the device and unplug the power cable from the power socket.
7. Make sure that the air tube contains no condensation. In case of condensation, switch on the device to dry up any moisture in the tube and then disconnect it as described above.
8. Empty out the remaining medicine from the nebulizer and clean the device in all of its parts, as described in the “Cleaning and Disinfection” section.

This device has been designed to operate 20 min. ON / 40 min. OFF.

Do not cover the device with a cloth, towel or anything else when in operation to avoid overheating or malfunction.

CLEANING AND DISINFECTION

Carefully clean all parts and remove any medicine residues and possible impurities after each treatment. The compressor (1) and the air tube (11) must be wiped clean with a clean damp cloth.

 Always wash hands thoroughly before cleaning and disinfecting the accessories.

 Wash/sterilize every part of the nebulizer.

 Do not expose the compressor to water or heat.

 Replace the air tube when passing to a new patient or in case of impurities.

 **Before cleaning, disconnect the plug from the electric outlet.**

 **The compressor is not waterproof and must not come into contact with liquids. In case of contact with liquids, immediately disconnect the power cable, shut off the device and dry it with care. Do not immerse the device in liquids.**

Cleaning with water

- Clean all of the device parts (**except for the air tube**) under hot running water for about 5 minutes, adding a small amount of mild detergent if necessary, following the dosage and limits specified by the detergent manufacturer.
- Rinse thoroughly, taking care to remove all residues, then leave to dry.

Disinfection

- All parts (**except for the air tube**) can be disinfected with chemical disinfectants, using doses and limits as specified by the manufacturer of the disinfectant.
- Disinfectants are generally available from pharmacies.

Steam sterilisation

- All of the components (**except for the air tube and the masks**) can be steam-sterilised up to a temperature of 100°C (15-20 min.).
- After sterilising, always leave the components to cool down to ambient temperature before using again.

 Do not repeat the sterilisation cycle when the components are still hot.

Sterilisation in boiling water

- All of the components (**except for the air tube and the masks**) can be immersed in boiling water for no longer than 1 minute.

 Boil the components in plenty of water.

 Disassemble all components of the nebulizer and boil them separately.

MAINTENANCE, STORAGE AND CARE

Order all spare parts from your retailer and trusted pharmacy or contact COLPHARMA customer service (see intro).

- It is recommended to replace the nebulizer, mouthpiece, nosepiece and masks after 1 year of use.
- Do not dry the nebulizer in the oven or with a hair-dryer or its components may be damaged.
- Always check that the filter is clean and replace it immediately if dirty or after a maximum of 60 days of use. Spare filters are provided with the device.
- To replace the filter (14), remove the filter cover from its seat. Remove the used air filter and insert the new filter. Reposition the filter cover in its seat.
- Do not use the device without a filter installed. Use only original filters.
- Do not wash the filter. If the filter is wet, replace it. If the filter is damp it may cause the device to stop operating.
- Wash the filter cover regularly (do not boil it) then dry it carefully before repositioning it in its seat.

 Never open the compressor compartment.

- Repairs must be performed exclusively by specialized personnel.

TROUBLESHOOTING

The device does not turn on

- Make sure that the plug is well inserted into the wall socket.
- Make sure to have pressed the ON/OFF (2) button.
- Make sure that the device is operating within the limits established in this instruction manual (20 min. ON / 40 min. OFF).

Nebulization is weak or non-existent

- Make sure that the air tube (11) is correctly connected at both ends.
- Make sure the air tube is not squashed, bent, soiled or blocked. If necessary, replace with a new one.
- Make sure the nebulizer is completely assembled (13) and ensure the components are correctly assembled.
- Make sure the medicine has been poured into the nebulizer.
- Make sure the filter is clean and the filter cover is positioned correctly.

The device is noisy

- Make sure to have inserted the filter correctly in its seat (14).

Do not attempt to repair the device. Contact COLPHARMA customer service.

WARRANTY

This device has a **2-year** warranty from the date of purchase. The warranty is valid only with the original purchase receipt.

- The 2-year warranty is valid only for the compressor.
- Spare parts, such as nebulizer, masks, mouthpiece, nosepiece, air tube and filters are guaranteed against defects in materials and workmanship for a period of one year. The warranty does not cover deterioration caused by normal wear and tear.
- Opening or tampering with the device voids the warranty.
- The warranty does not cover damage caused by improper use, accidents or failure to abide by the instructions for use.

TECHNICAL SPECIFICATIONS

Nebulizer:	EASY 2
Maximum filling volume:	8 ml
Minimum filling volume:	2 ml
Residual volume:	0.5 ml
Nebulization speed (NaCl 0.9%):	0.35 ml/min
Particle size:	4.07 μm – 61% < 5 μm


Power supply:	AC 230 V/50 Hz, 190 VA
Maximum pressure:	1.35 bar
Operational flow:	5 L/min
Sound level:	55 dBA
Operating cycle:	20 minutes ON, 40 minutes OFF
Operating conditions:	+5 °C – +40 °C; 15% – 90% R.H.; 700 – 1060hPa
Shipping and storage conditions:	-25 °C – +70 °C; 15% – 90% R.H.; 700 – 1060hPa

Product weight:	1.1 kg
Dimensions:	18 × 13.2 × 9 cm
Pollution degrees:	Degree 2
Altitude:	High altitude up to 2000 m
Useful life:	1000 hours
Power supply:	mains plug
IP class:	IP21

We reserve the right to apply technical changes.

Particle dimensions diagram

This diagram might not be applicable in the case of suspensions or highly viscous medicines. More information is available from the manufacturer of the medicine in question.

MANUFACTURER'S GUIDE AND DECLARATION - ELECTROMAGNETIC EMISSIONS - FOR ALL EQUIPMENT AND SYSTEMS

Warning!

- The use of this equipment adjacent to or stacked with other devices is to be avoided as it could adversely affect correct operation. If this type of use is necessary, this equipment and other devices must be monitored to ensure that they are all operating correctly.
- The use of accessories, transducers or cables other than those specified or supplied by the manufacture of this appliance may lead to an increase in electromagnetic emissions or a reduction in electromagnetic immunity of the appliance, causing improper operation.
- Portable RF communication equipment (including peripherals such as cables for antennas and external antennas) must be used at least 30 cm (12 inches) from any part of the equipment, including the cables specified by the manufacturer. If this is not the case, it may lead to deterioration in equipment performance.
- In the test conditions specified for immunity, the product is able to guarantee basic safety and essential performance.
- The emission characteristics of this equipment make it suitable for use in home healthcare (CIS-PR11 Class B).
- The switch serves to cut off the device from the mains electricity; its direction complies with IEC 60447.

INFORMATION ON CONFORMITY FOR ALL EMC TESTS

Electromagnetic emissions (home healthcare environment)

Emissions test (IEC60601-1-2:2014)	Conformity
Radiated and conducted emissions (RF emissions)	CISPR 11 group 1 class B
Harmonic emissions (IEC 61000-3-2)	Class A
Voltage fluctuations/flicker (IEC 61000-3-3)	Compliant

Declaration - Electromagnetic immunity (home healthcare environment)

Immunity Test	Test Level	Conformity Level
Conducted RF (IEC 61000-4-6:2013)	3V 150kHz at 80 MHz 6 V in ISM and between 0.15 MHz and 80 MHz	3V 150kHz at 80 MHz 6 V in ISM and between 0.15 MHz and 80 MHz
Radiated RF (IEC 61000-4-3:2006+A1:2007+A2:2010)	10V/m 80 MHz at 2.7 GHz	10V/m

Declaration - Electromagnetic immunity (home healthcare environment)

Immunity Test	IEC 60601 test level	Conformity level
Electrostatic discharge (ESD) IEC 61000-4-2:2008	±8kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transients/bursts IEC 61000-4-4:2012	±2kV per power supply line	±2kV per power supply line
Overvoltage IEC 61000-4-5:2005	±0.5kV, ±1 kV line to lines	±0.5kV, ±1 kV line to lines
Voltage drops, brief power outages, and voltage variations on the electrical supply lines IEC 61000-4-11:2004	0% UT, 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT, 1 cycle and 70% UT, 25/30 cycles Single phase at 0° 0% UT, 250/300 cycles	0% UT, 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT, 1 cycle and 70% UT, 25/30 cycles Single phase at 0° 0% UT, 250/300 cycles
Magnetic field and mains power frequency (50/60Hz) IEC 61000-4-8:2009	30 A/m	30 A/m

NOTE: The EUT is the mains voltage in AC before application of the test level. The following phenomenon meets basic safety and essential performance requirements.

UT: 230V~/50Hz – The EUT shifts from the normal value but remains above 10 psi when the flow is 3 l/min.

UT: 230V~/50Hz – The EUT ceases to function when 0% UT (total absence of voltage) is applied, but can automatically reset to normal mode automatically.

Declaration – immunity to proximity effects from RF wireless communication equipment

Immunity test	Test frequency	Modulation	Maximum power	Immunity level	Conformity level
Radiated radio frequency IEC 61000-4-3: 2006+A1:2007 +A2:2010	385 MHz	**pulse modulation: 18 Hz	1.8W	27V/m	27V/m
	450 MHz	*FM + 5Hz deviation: sinusoid at 1kHz	2W	28V/m	28V/m
	710 MHz	**pulse modulation: 217 Hz	0.2W	9V/m	9V/m
	745 MHz				
	780 MHz				
	810 MHz				
	870 MHz	**pulse modulation: 18 Hz	2W	27V/m	27V/m
	930 MHz				
	1720 MHz				
	1845 MHz				
	1970 MHz	**pulse modulation: 217 Hz	2W	28V/m	28V/m
	2450 MHz				
	5240 MHz	**pulse modulation: 217 Hz	0.2W	28V/m	28V/m
	5500 MHz				
	5785 MHz				

NOTE * – As an alternative to FM modulation, it is possible to use a pulse modulation of 50% at 18 Hz, since although it does not represent real modulation, it would constitute a worst-case scenario.

NOTE ** – The carrier must be modulated using a square wave signal with 50% work cycle.

MATERIAL COMPATIBILITY

There is no specific or unusual risk associated with the compatibility between the materials of the component and the liquid dispensed.

The materials used in the nebulizer kit might not be compatible with solutions/suspensions/emulsions that have not been assessed.

The phthalates conform to European REACH regulations.

INFORMATION ON DISPOSING OF THE DEVICE

Dispose of the device, components and optional accessories in compliance with applicable local regulations.

Unauthorised disposal may cause environmental pollution.

REGULAR SAFETY CHECKS

Inspection and preventive maintenance to perform, including the frequency of this maintenance.

1. Clean the power cable plug at least once a year. An excessive build-up of dust on the plug is a fire hazard.
2. The following safety checks must be performed at least every 24 months by a technical engineer from the manufacturer who is suitably trained, and with the knowledge and practical experience needed to perform these tests.

- a) Inspect the equipment and accessories for any mechanical or operational damage.
- b) Make sure that all safety labels are clearly legible.
- c) Inspect the fuse to ensure its conformity with the rated current and breaking capacity.
- d) Make sure that the device operates correctly as described in the instructions for use.
- e) Test the leakage current of the casing according to IEC 60601-1.

Limit NC 100 μ A, SFC: 500 μ A.

- f) Test the leakage current of the patient according to IEC 60601-1.

Limit: for alternating current (AC): 100 μ A, for direct current (DC): 10 μ A.

- g) Test the leakage current of the patient in single fault condition according to IEC 60601-1.

Limit: for alternating current (AC): 0.5 mA, for direct current (DC): 50 μ A.

Leakage current must never exceed the limit. The figures must be recorded in a register for the appliance.

If the device does not operate correctly or does not pass one of the above tests, it must be repaired.

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