COLPHARMA

Air 2000

P R O

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EN	1. INSTRUCTIONS	page 22
	2. SAFETY INFORMATION	page 23
	3. PRECAUTIONS	page 24
	4. PREPARING FOR USE	page 24
	5. ASSEMBLING THE DEVICE	page 25
	6. USING THE DEVICE	page 26
	7. CLEANING AND DISINFECTION	page 27
	8. MAINTENANCE AND STORAGE	page 29
	9. NOTES	page 30
	10. DECLARATION	page 30
	11. TROUBLESHOOTING	page 31
	12. WARRANTY	page 32
	13. CLASSIFICATION	page 32
	14. MANUFACTURER'S DECLARATION	page 33
	15. CONFORMITY INFORMATION FOR EACH ELECTROMAGNETIC COMPATIBILITY TEST	page 34
	16. TECHNICAL SPECIFICATIONS	page 35
	17. REGULAR SAFETY CHECKS	page 36
	18. DISPOSAL	page 36

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Portable nebuliser with EASY MESH technology

- | | |
|-----------------------------|---|
| 1. Nebuliser unit | 5. Nebulisation speed adjustment button |
| 1a. Lid | 6. Operating LED |
| 1b. Fastening tab | 7. Battery compartment cover |
| 1c. Medication container | 8. C-type USB port |
| 1d. Replaceable Easy Mesh | 9. Mouthpiece |
| 1e. Fitting for accessories | 10. Adult mask |
| 1f. Fastening clip | 11. Child mask |
| 2. Power contacts | 12. Type-C USB port |
| 3. Power unit | 13. Batteries (2 x AA) |
| 4. ON/OFF button | |

The patient needs to choose the most suitable mask size to suit their face.

The mouthpiece (9) and masks (10 and 11) are not disposable accessories. Since these items can deteriorate with increasing uses, it is recommended that you replace them every six months.

1. INSTRUCTIONS

Dear Customer,

Thank you for choosing to purchase an AIR 2000 PRO from Colpharma, the portable nebuliser with EASY MESH technology (replaceable mesh).

Intended use: this product is a mesh nebuliser treatment system, designed to turn medications into a fine aerosol to be inhaled by the patient.

Recommended use: this device has been developed for the effective treatment of asthma, allergies and other respiratory problems.

Operating principle: when the device is switched on, the high-frequency vibrations generated by a circuit cause the ceramic plate and mesh to oscillate at high frequency forcing the medicine to flow through the microscopic holes in the mesh and be nebulised into ultra-fine particles to be inhaled by the patient.

Intended users: the device can be used on both adults and children.

Usage environment: the device is intended for use in home, clinic, and hospital environments.

Contraindications: none.

2. SAFETY INFORMATION

READ ALL THE INSTRUCTIONS CAREFULLY BEFORE USING THIS DEVICE.

Like all medical devices, this product can become unusable in the event of a power outage, if the batteries should become flat, or in the event of a mechanical problem. We recommend keeping a spare device handy.

When using electrical products, always follow the basic safety precautions. Like for all electrical devices, take particular care around children.

WARNING

- Only use this device with medicines prescribed by a physician and follow the dosage, duration and frequency instructions for the prescribed treatment.
- This nebuliser is exclusively for use in respiratory therapy; any other use of the device is improper and therefore, dangerous. The manufacturer cannot be held responsible for any other damage caused by improper or incorrect use.
- Do not share the nebuliser with other people. This product (mask and mouthpiece included) must be used by a single patient.
- Clean and disinfect the nebuliser unit and accessories (mouthpiece and mask) as described in the section "Cleaning and disinfection". Make sure that the nebuliser unit and accessories have been cleaned and disinfected correctly before use, to avoid all possible contamination.
- The patient is the intended operator of the nebuliser. Adult supervision is necessary if the device is being used by a child or person requiring special assistance.
- For safety reasons, always disconnect the Type-C USB cable from the Type-C USB port in the following situation:
 - if the power adaptor is damaged
 - if a malfunction has occurred during operation
 - before cleaning the device
 - immediately after use.
- The device must not be used if there is any suspicion of a breakdown.
- Do not store the device and its accessories in a damp area. Contamination and residual moisture can favour the spread of bacteria, increasing the risk of infection.
- Never handle the power adapter with wet hands.
- Never use the device if the Type-C USB cable or power adaptor are wet.
- Never remove the power adaptor from the main socket by pulling the cable.
- Do not swallow small parts of the nebuliser.
- This device is for human use only.
- To avoid all risk of entrapment and strangulation, keep the cables out of the reach of small children.

3. PRECAUTIONS

- Do not try to clean the mesh nebuliser with foreign objects, as there is a risk of damage.
- Never use this device when the environmental temperature is higher than 40 °C. For more requirements pertaining to the operating environment, see the "Technical Specifications" section.
- Keep the device out of direct sunlight, and away from excessive heat or cold, to prevent damage to the batteries.
- Do not shake the nebuliser during operation.
- Avoid subjecting the power unit and its component parts to high impact, for example, such as falling to the ground.
- Never try to open, repair or modify this device.
- Only use the accessories named in this instruction manual.
- Before use, make sure that the device and all its parts and components (for example, mouthpiece and mask) are correctly assembled, as described in these instructions for use.
- Comply with local laws and recycling programmes regarding the disposal or recycling of components, batteries and packaging.
- Never use water to wash the power unit, adaptor or Type-C USB cable. If any liquid is spilled on the power unit or adaptor, clean them immediately.
- Never place the power unit in water or other liquids, and be sure to protect it from splashes or sprays of water or other liquids.
- Do not use the device in damp environments such as bathrooms.
- When connecting the power unit, make sure that the voltage at the power socket is correct and that the plug does not cause an overload.
- Disconnect the power after using the device, and never leave a device unattended if under tension.
- If the ambient temperature where the device is stored does not fall within an interval of 5-40 °C, leave it to stand for at least 30 minutes before use.

4. PREPARING FOR USE

POWER SUPPLY

Meaning of the different colours of operating LED:

- 1) Whether the device is powered by battery or Type-C USB cable, if the operating LED lights up green, this means that the device is ready for use.
- 2) In battery mode, when the battery charge is too low to guarantee the normal operation of the device, the operating LED (battery indicator flat), will turn orange without blinking.
- 3) When the operating LED flashes green, the device enters cleaning mode.

PREPARING THE BATTERIES OR C-TYPE USB CABLE

1) Using with batteries

Open the lid to the battery compartment (bottom of the power unit) following the direction of the arrows (Fig. A). Fit two AA batteries in the battery compartment and close the lid.

Note: Follow the (+) and (-) signs on the lid of the battery compartment to refit the batteries in the correct direction.

Battery duration and replacement

- 1) The nebuliser can operate using alkaline or rechargeable nickel-metal hydride (NiMH) batteries
- 2) When the operating LED turns orange, replace both batteries with new ones.
- 3) The duration of the batteries depends on their capacity and their condition.

2) Use with Type-C USB cable

- 1) The Type-C USB cable has two ends: one is type C and the other is USB. Connect the Type-C end to the Type-C port on the right side of the power unit (see Fig. B).
- 2) Connect the other USB cable end to the power unit. The power unit must be able to provide a direct current of 5V/1A and may be a medical adaptor or power bank.

OPERATING TEST

If using batteries, switch on the device and make sure that the operating LED switches on.

If this is not the case, check the batteries and replace them if necessary.

If using the Type-C USB cable, switch on the device and make sure that the operating LED switches on. If this is not the case, check the power supply unit and replace it with a suitable unit that meets all requirements.

Notes:

- NEVER insert a battery with the positive pole (+) and negative pole (-) facing the wrong way.
- NEVER use old and new batteries together or mix and use different battery brands together.

Select a medical power adaptor that complies with the regulation IEC/EN 60601-1.

5. ASSEMBLING THE DEVICE

Fitting the nebuliser unit:

1. Push the nebuliser unit in the direction shown by the arrow (Fig. C) in the slot of the power unit. A "click" will indicate that the nebuliser is correctly in position. Engage the fastening clip (1f) vertically in the direction shown in the figure, ensuring that the two cylindrical projections are inserted into the corresponding circular holes on the fitting for accessories (1e) (see Fig. D).
2. If the device is not being used immediately, it must be stored in its bag, clean, dry and assembled.

6. USING THE DEVICE

1. Opening the medication compartment

To open the medication compartment, hold the power unit in one hand and grip the tab on the lid between the thumb and forefinger of the other hand to release it by pulling upwards (Fig. E).

2. Adding the medication

Make sure that the power LED is off before adding the medication. Add the amount of medication prescribed by the physician (Fig. F).

The maximum and minimum volume of the fluid is 8 ml (to the MAX LEVEL mark) and 0.5 ml respectively. Never exceed the LIMIT (max. 8 ml).

3. Closing the medication container

Close the lid of the medication container (Fig. G). The lid should click into place if closed correctly.

4. Fitting the accessories and proceeding to inhalation

Fitting the mouthpiece: connect the mouthpiece to the connection of the nebuliser unit (Fig. H).

Using the mouthpiece: bring the mouthpiece up to the mouth and gently seal the lips around it to allow the nebulised medication to flow freely (Fig. I). Press the ON/OFF button.

Fitting the mask: connect the mask to the fitting on the nebuliser, following the direction of the arrow (Fig. J).

Using the mask: bring the mask up to the face (Fig. K) and adjust the elastics to secure it comfortably. Press the ON/OFF button.

- 1) In case of administering medication through the nose, it is recommended to breathe in deeply through the nose and to breathe out through the mouth.
- 2) In case of administering medication via the throat, it is recommended to breathe in deeply through the mouth and to breathe out through the nose.
- 3) It is possible to alternate the two methods described above. If using the mask, gently press it onto the face so that it perfectly adheres to nose and mouth.

Notes:

- For reasons of safety, always use original accessories to replace old ones.
- Use the mouthpiece correctly to prevent a build-up of medication on the tongue.
- Do not swallow the nebulised medication and make sure it reaches the affected area.
- For effective aerosol treatment, it is necessary to be calm, relaxed and seated. An uncomfortable posture and uneven breathing may cause obstruction in the airways.
- Be sure to be seated in the correct position; breathe deeply, evenly, and slowly to help the nebulised medication to be inhaled into the airways.
- It is possible to press the nebulisation speed adjustment button to change the nebulisation flow. The device has three nebulisation speed levels; the initial level is maximum.
- Press the ON/OFF button to interrupt the aerosol therapy, when required. The device switches off automatically 20 minutes after being switched on.
- When there are lots of bubbles between the medication and the mesh, the device will not be able to nebulise the medicine. It is possible first to switch off the device then to shake it gently.

Then restart the device (Fig. L).

- To guarantee the normal nebulisation of the medication, keep the device upright or tilt it slightly towards the mesh.
- Slightly tilt the nebuliser towards you when the solution is almost finished to make sure that any remainder stays in contact with the mesh to be nebulised completely.
- A high-viscose solution may not nebulise fully or may even clog the mesh. In this case, switch off the device and use distilled water to rinse out any built up solution.
- During use, take care not to shake the device since this could cause malfunctions, such as automatic switch off.

5. End of inhalation

At the end of the treatment, press the ON/OFF button to switch off the device. Disconnect the Type-C USB cable from the power unit and disconnect the adaptor from the socket if the Type-C USB cable was used as a power supply.

Remove the mouthpiece or mask from the device, following the assembly steps in reverse.

Pour away any residual medication and clean the device. For more information, see the "Cleaning and Disinfection" section.

7. CLEANING AND DISINFECTION

WARNING

The nebuliser unit and accessories (mouthpiece or mask) must be cleaned and disinfected after every use to prevent the spread of micro organisms which increase the risk of infection.

LIMITS TO USE

Accessories can be used up to 300 times.

PREPARATION

Remove the mouthpiece or mask from the nebuliser unit after every use.

CLEANING

Clean the nebuliser unit after every use.

1. Remove the mouthpiece or mask from the nebuliser unit after every use.
2. Open the lid on the nebuliser unit to eliminate any residual medication.

Add tap water to the medication container and gently shake it, then pour out the water. Take care not to wet the outer surface of the device.

3. **Activating cleaning mode:**

Add more tap water and close the nebuliser lid. Press and hold down the ON/OFF button for 5 seconds; the LED will flash green and the device will enter the cleaning mode. After 2-3 minutes, press the ON/OFF button to interrupt the cleaning mode (if no operation is performed during the cleaning process, the device will switch off automatically after 10 minutes). Pour out the remaining

water at the end of the process. During this process, the cleaning mode will be interrupted in case of low battery charge.

4. Remove the nebuliser unit from the power unit. Disconnect the Type-C USB cable from the power unit if using the USB power supply. Wash and rinse the nebuliser unit with distilled water (the power unit must not be washed). Shake out the excess water and leave to dry naturally on a clean, dry towel.

Do not touch the mesh with fingers, cotton buds or other items.

5. If necessary, use a piece of gauze or a clean towel to remove any marks from the power unit.
6. Make sure that all clean parts are completely dry before putting them away or using them.

Warning:

- The power unit contains electronic parts and electrodes. Do not immerse or rinse the power unit in water. There is a risk of electric shock or loss of current, with resulting malfunction of the device.
- Do not touch the mesh with fingers, cotton buds or any other items.
- Do not wash device parts in the dishwasher.
- Never use a microwave to dry parts.
- When the mesh is seriously clogged, it cannot be cleaned with the above methods; it is possible to try the following method: add boiling water and a few drops of white vinegar to the medication container and activate the cleaning mode. If the mesh remains clogged, replace the part.

HOW TO REPLACE THE MESH (EASY MESH TECHNOLOGY)

1. Detach the nebuliser unit from the power unit by pressing the PRESS button with one hand and sliding the nebuliser unit in the direction indicated by the arrow with the other hand (Fig. M).
2. Take out the fastening clip, vertically and turn the fitting for accessories anticlockwise in the direction shown by the arrow, (Fig. N).
3. Press your finger into the two raised circular holes on the back of the nebuliser unit to push the Mesh outwards towards the fitting for accessories. This will remove the Mesh so that it can be replaced (Fig. O).
4. Insert a new Mesh and refit the fitting for accessories by following the steps in reverse order to those described above (Fig. P), ensuring that the fastening clip is fitted at the end.

CLEANING MASKS AND MOUTHPIECE

1. Immerse all the disassembled parts in warm tap water (40-50 °C) for at least 10 minutes. If necessary, use a clean brush to remove excess dirt (the brush must be reserved exclusively for this purpose).
2. Carefully rinse all the parts in running water (the flow rate must be above 3 litres per minute) for at least one minute each.
3. Independently of the whether or not the parts are disinfected immediately, dry all pieces with a fresh piece of medical gauze. It is possible to remove the excess water rapidly by shaking the parts.

**Warning:**

Insufficient disinfection helps to spread bacteria and therefore increases the risk of infection. Suitable disinfection can only be guaranteed if the application times are respected and if all individual parts are completely immersed in the solution for the full application time. There must not be any cavities or air bubbles.

DISINFECTION

After cleaning, disinfect the disassembled nebuliser unit and the accessories used (mouthpiece or mask). Efficient disinfection is only possible if the nebuliser unit and accessories have been cleaned. Use a solution of hydrogen peroxide at 2.0% (p/p) to disinfect the nebuliser unit and the accessories. To guarantee safety while handling chemical products, follow the user instructions for the disinfectant and in particular, the accompanying safety instructions.

8. MAINTENANCE AND STORAGE

- Make sure there are no liquid residues in the nebuliser unit.
- Make sure that the nebuliser unit is cleaned and disinfected carefully after every use.
- Store the device and accessories in a clean, safe place.
- Protect the device from direct sunlight or knocks.
- Remember to keep any heating device or naked flames well away from the storage place.
- Protect the device from contact with corrosive fluids.
- Never use a microwave oven to dry the device or its accessories.
- Never wind the power cord around the device.
- If repairs are necessary, always contact your supplier or nearest local retailer.
- This device is provided with a five-year warranty, valid from the date of purchase.
- Only use or store the device in an environment where the temperature and moisture levels are compliant with the specific requirements in the "Technical Specifications", otherwise the device may not operate correctly.

**Warning:**

Do not disassemble the device or attempt to repair it without authorisation before returning it to the manufacturer for inspection and repair.

9. NOTES

- This device is not suitable for use in anaesthesia breathing systems or in pulmonary ventilation systems. This device does not require a gas supply.
- Never touch the metal part of the mesh in the nebuliser unit with your hands.
- Never use fluids containing fat or suspended particles (including herbal juices for diffusion); please note that specifications may vary according to the type of medication used; please use the standard fluids recommended by a physician for nebulisation.
- Nebuliser solution must be prepared using a 0.9% saline solution or as recommended by the physician. Tap water or boiled water must not be used.
- The container to prepare the medication must be sterilised by boiling.
- The device must not be used for longer than 15 minutes at a time or than any other limits indicated by the physician.
- Considering that there is a risk that children or babies could swallow small parts, the device may only be used on children or babies with the appropriate supervision. Children and patients affected by mental illness can only use this device with the supervision of a responsible adult.
- Remove the batteries if the device is not to be used for an extended period.
- Remove the batteries if the device is going to be taken away from the house or when travelling.
- Do not place the device in water, under running water or in very damp environments.
- Never use the device in the shower.
- If the device is not working, see the "Troubleshooting" section to find the solution.
- Always use personal accessories, as stated in this document.
- Never use the device in places where it could be exposed to flammable gases.
- This device conforms to the EMC IEC 60601-1-2 standard. Nevertheless, if used together with other medical devices or electrical equipment, these could affect the operation of one of the devices; please do not use this device with other medical devices or electrical equipment.

10. DECLARATION

If the technical information concerning the user's needs do not fall within the Company's confidentiality, this latter undertakes to provide this information in compliance with the envisaged procedure; this includes wiring diagrams, lists of parts, and other related technical information that is not covered by trade secrets. As far as the channels and procedures for accessing information are concerned, please contact your retailer or the manufacturer.

11. TROUBLESHOOTING

No.	Problem	Possible causes	Solution
1	Insufficient nebuliser speed	The battery flat indicator switches on / flat battery.	Replace the batteries and restart the device.
		The mesh holes are clogged.	Clean the nebuliser (check the "Cleaning and Disinfection" section). Replace the mesh, if necessary.
		The medication is finished or out of contact with the mesh for more than 10 seconds.	Tilt the device towards the user to allow the medication to come into contact with the mesh.
		The mesh is damaged.	Replace the mesh.
2	The operating LED does not switch on and the device is not working	The batteries are inserted facing the wrong way.	Follow the (+) and (-) on the lid of the battery compartment to refit the batteries in the correct direction.
		Low battery charge.	Replace the batteries or use the Type-C USB cable to connect to the power supply or power unit.
		The C-Type-C USB cable is not correctly connected to the power unit or supply.	Replace the power supply or Type-C USB cable.
		Lack of power or C-Type-C USB cable damaged.	Refit the nebuliser unit correctly.
3	The operating LED switches on but the device is not working	Incorrect nebuliser unit installation. The battery flat indicator switches on / flat battery.	Refit the nebuliser unit correctly.
		The power contacts on the power unit or those in the nebuliser unit are dirty.	Replace the batteries and restart the device.
		The angle of the device is too great and the medication cannot come into contact with the mesh.	Clean the power contacts with a medical gauze.
			Reduce the tilt angle of the device or tilt it towards the user to allow the medication to come into contact with the mesh.

4	The nebuliser switches off automatically during use	The angle of the device is too great and the medication cannot come into contact with the mesh.	Reduce the tilt angle of the nebuliser or tilt the front part towards the user to allow the medication to come into contact with the mesh.
		The nebuliser is being shaken roughly during use.	Hold the nebuliser with care during use.

If failures persist after trying any of the measures above, do not try to repair the device; contact your retailer or the COLPHARMA customer sales service.

12. WARRANTY

This device has a **five-year warranty**, valid from the date of purchase. The warranty is valid only with the original purchase receipt.

- The warranty covers the device against defects in materials and workmanship.
- Consumable parts, such as the nebuliser unit and replaceable mesh are covered by a warranty against manufacturing defects upon first use, but not thereafter, given the natural wear and tear that occurs over time.
- The spare parts such as masks or mouthpiece, are guaranteed for one year against defects in materials and workmanship, but not against deterioration caused by normal wear and tear.
- Opening or tampering with the device voids the warranty.
- The warranty does not cover damage caused by improper use, accidents or failure to abide by the instructions for use.

WARNING

- Never make changes to the device without the manufacturer's authorisation.
- All important maintenance interventions on the device must be performed by an authorised assistance service or retailer. No parts inside the device can be repaired by the user; before any intervention, contact an authorised dealer or the manufacturer.

13. CLASSIFICATION

Appliance not suitable for use in the presence of flammable mixes.

- BF type part applied. The mask or mouthpiece are recognised as an applied part.
- Appliance powered internally (without adaptor), Class 1 equipment (with adaptor).

14. MANUFACTURER'S DECLARATION

The software development process for NB-818B conforms with the EN 62304 standard.

Directives and manufacturer's declaration - electromagnetic immunity

The mesh nebuliser, model NB-818B is intended for use in the electromagnetic environment specified below. The client or user of the mesh nebuliser model NB-818B must ensure that it is used in this environment.

Immunity test	Test level IEC 60601	Conformity level
Electrostatic discharge (ESD) IEC 61000-4-2:2008	±8kV in contact ±2 kV, ±4 kV, ±8 kV, ±15 kV in air	±8kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV in air
Electrical fast transients/burst IEC 61000-4-4:2012	±2kV per power supply line	±2kV per power supply line
Overvoltage IEC 61000-4-5:2005	0.5kV, ±1kV, ±1 kV line to lines	±0.5kV, ±1 kV line to lines
Voltage dips, short interruptions and voltage variations in the power network IEC 61000-4-11:2004	0% UT, 0.5 Cycles at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT, 1 Cycle and 70% UT, 25/30 cycles Single phase: at 0° 0% UT, 250/300 cycles	0% UT, 0.5 Cycles at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT, 1 Cycle and 70% UT, 25 cycles Single phase: at 0° 0% UT, 250 cycles
Supply frequency (50/60 Hz) magnetic field IEC 61000-4-8:2009	30 A/m	30 A/m
Conducted RF (IEC 61000-4-6:2013)	3 V 150 kHz at 80 MHz 6 V in the ISM band and in the band between 0.15 MHz and	3 V 150 kHz at 80 MHz 6 V in the ISM band and in the band between 0.15
Radiated RF (IEC 61000-4-3: 2006+A1: 2007+A2:2010	80 MHz 10 V/m 80 MHz at 2.7 GHz	MHz and 80 MHz 10V/m

Portable and mobile radio-frequency communication equipment (including peripherals such as aerial cables and external aerials) cannot be used at a distance of less than 30 cm from any part of the device, including the cables specified by the manufacturer. If this is not the case, it may lead to deterioration in equipment performance.

NOTE 1: at 80 MHz and 800 MHz, the higher frequency range is used.

NOTE 2: these guidelines may not be applicable in all situations. The electromagnetic propagation is influenced by power draw or by reflection by structures, objects or persons.

- a. The field intensity generated by fixed transmitters such as for RT and MT base stations, radio amateurs, AM and FM radio transmissions, TV transmissions, cannot be precisely envisaged at a theoretical level. To assess the electromagnetic environment caused by fixed RF transmitters, it is appropriate to consider an electromagnetic survey in situ. If field intensity measured in the area in which the mesh nebuliser, model NB-818B exceeds the applicable RF conformity level as above, it is necessary to observe the mesh nebuliser, model NB-828B to ensure normal operation.
If any abnormal performance is found, additional measures may be required, such as the reorientation or positioning of the mesh nebuliser model NB-818B.
- b. In the frequency range from 150 kHz to 80 MHz, field intensity must be less than 3 V/m.

15. CONFORMITY INFORMATION FOR EACH ELECTROMAGNETIC COMPATIBILITY TEST

Electromagnetic emissions (home healthcare environment)	
Overvoltage IEC 60601-1-2:2014	Conformity
Radiated and conducted RF emissions	CISPR 11 group 1 class B
Harmonic emissions (IEC 61000-3-2)	Class A
Voltage fluctuations/flicker (IEC 61000-3-3)	Compliant

- Use of this equipment is to be avoided close to or stacked on other equipment, since this could lead to incorrect operation. If such use is required, make sure that this and the other equipment are operating correctly.
- The use of accessories, transducers or cables other than those specified or supplied by the manufacture of this appliance may cause an increase in electromagnetic emissions or a reduction in electromagnetic immunity of the appliance, causing incorrect operation.


16. TECHNICAL SPECIFICATIONS

Product name:	Portable mesh nebuliser
Model:	NB-818B
Power supply:	Adaptor (CA-CC 5 V, 1 A) or 2 AA batteries (3V)
Rated power:	<3 W
Ultrasonic frequency:	Approx. 110 kHz
Safety level:	Type BF
Maximum nebuliser capacity:	≥0,2 ml/min
Maximum/minimum fluid volume:	8 ml/0,5 ml
MMAD:	≤ 5 µm
FPD:	> 60%
Flat battery indicator:	1.9V±0.2V
Duration:	Main unit: 5 years
	Mesh: 180 hours
	Accessories: 6 months
Battery duration:	90 minutes with two AA alkaline batteries (or LR6)
Self-cleaning mode:	Yes
Operating temperature and humidity – Air pressure:	+5 °C - +40 °C, 15% - 90% relative humidity, from 700 hPa to 1060 hPa
Storage temperature and humidity - Air pressure:	-25 °C - +70 °C, 15% - 90% relative humidity, from 700 hPa to 1060 hPa
Pollution level:	Grade 2
Voltage category:	Category II
Altitude:	<3000 m
Size:	108 mm × 75 mm × 45 mm

Technical specifications may undergo variations without warning. Performances may vary with medicines such as suspensions or high viscosity. For more details, see the technical data sheets for the medicine suppliers.

Pursuant to regulation EN 13544-1: 2007+A1:2009 "Respiratory therapy equipment - Part 1: Nebuliser systems and relevant parts", annex CC.3, which includes the use of a multi-stage cascade impactor to measure particle dimensions.

Test conditions	
Test temperature	24 ±2°C
Test humidity	45%-75% R.H.
Test solution	NaF 2,5% (M/V)

17. REGULAR SAFETY CHECKS

If using the device with an adaptor, it is necessary to perform preventive inspections and maintenance interventions, taking the frequency of these interventions into account.

- Before each use, check the adaptor; do not use if damaged.
- Clean the adaptor plug at least once a year. An excessive build-up of dust on the plug is a fire hazard.

In the interest of progress, the manufacturer reserves the right to make technical changes without prior warning. No prior warning will be given in the event of changes made to this manual. The trademarks and names mentioned are the property of their respective owners.

18. DISPOSAL

For environmental reasons, never dispose of the device in household waste. Please dispose of the appliance according to the EC WEEE Directive (Waste electrical and electronic equipment). In case of any questions, please contact your local authority for information on waste disposal.

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HONSUN (NANTONG) Co., Ltd
Shanghai International Holding Corp. GmbH (Europe)