Version:Ver.1.4

ENSM80-8420306019WHH /000000000000

Applicable Models: YK-80A, YK-80B, YK-80C, YK-81A, YK-81B,

YK-81C, YK-82A,YK-82B, YK-82C, YK-83A, YK-83B, YK-83C,

YK-84A, YK-84B,YK-84C, YK-81CEU.

Xuzhou Yongkang Electronic Science Technology Co., Ltd.
16#2nd Floor,60#1,60#2,No.6 Building 1st Phase Economic Development
Manufacturing Zone,LIANGZI Valley, No.6 Leiyue Road Xuzhou
FUDU-231008 Xuzhou,PEOPLES REPUBLIC OF CHINA

Prolin GmbH
Brehmestraße 56, 40239 Düsseldorf
Germany
Email: med@eulin.eu
office@yonker.cn

	wavelength	radiation power
RED	Approx. 660nm	1.8mW
IR	Approx. 940nm	2.0mW

This parameter can be especially useful to clinicians.
5. PI index: measurement scope: 0%-20% (optional)
6. Power: two AAA 1.5V alkaline batteries
7. Automatic standby: the product shuts off by itself when no finger is in the sensor for about 8 seconds
8. Dimension: 58mm x 36mm x 33mm
9. Operation environment:temperature: 5℃ ~ 40℃
Humidity: 15%RH ~ 80%RH
Atmospheric pressure: 70kPa ~ 106kPa
Transport, storage environment:Temperature: -10℃ ~ 40℃
Humidity: 10%RH~95%RH
Atmospheric pressure: 70kPa ~ 106kPa, non-corrosive gas and well-ventilated environment.
10. Declaration: EMC of this product comply with IEC60601-1-2 standard.

11. Automatic start-up: The device has the function of SS automatic signal detection. When the finger is inserted, it will start up automatically.(optional)
12. Gravity sensing function: finger movement, the screen display will change with the gravity sensing changes (optional)
13. Service life: 5 years

Classification

1. Management Class for Medical Devices: II equipment
2.Anti-electric Shock Type: Internally powered equipment
3.Anti-electric Shock Degree: Type BF equipment

Maintenance and Preservation

1. Replace the batteries timely when low voltage lamp is on.
2. Clean the surface of fingertip pulse oximeter before it is used to diagnose patients.
3. Remove the batteries inside if you will not operate the Oximeter for a long time.
4. It would be better to preserve the product in -10~40℃ (-14~104 ℉) and humidity is 10%-95%.
5. It is recommended that the product should be kept dry anytime. A wet ambience might affect its lifetime and even damage the product.
6. Please follow the law of the local government to deal with used batteries.

Symbols and Definitions

	Type BF applied part		This device has no internal battery
	This way up		Humidity limitation

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Ver.1.4

Fingertip Pulse Oximeter

This Fingertip Pulse Oximeter is a kind of innovated medical device with non-invasive features for artery SPO2 and PR detection. Being portable, it is able to measure SPO2 and PR values quickly and precisely.

General Description

Hemoglobin Saturation is the percentage between the capacity of Oxyhemoglobin (HbO2) that compounded with oxygen and that of all combinative hemoglobin (Hb) in blood. In other words, it is the saturation of Oxyhemoglobin in blood. It is a very important physiological parameter for Respiratory and Circulation Systems. Many respiratory diseases could reduce hemoglobin saturation in human blood. Moreover, factors such as Automatic Organic Regulation Malfunction caused by anaesthesia, trauma resulted from major operation and some medical examination can also cause problems in oxygen supply, which might reduce human hemoglobin saturation. As a result, such symptoms as mesgrim, vomiting and asthenia might appear to patients. Hence, it is very important to know hemoglobin saturation of patient timely in clinical medical aspects.

The fingertip pulse oximeter features in small volume, low power consumption, convenient operation and portability. It is only necessary for patient to put one finger into fingertip photoelectric sensor for diagnosis, and the display screen will directly show measured value of hemoglobin saturation. It has been proved in clinical experiments that it possesses rather high precision and repeatability.

Measurement principle

The principle of the oximeter is as follows: An experience formula of data process is established by exerting Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive hemoglobinand oxyhemoglobin (HbO2) in glow and near-infrared zones. Operation principle of the instrument is to combine Photoelectric Oxyhemoglobin Inspection Technology with Capacity Pulse Scanning and Recording Technology, so that the different wavelengths (660nm glow and 940nm near infrared light) can be focused onto human nail through perspective clamp finger-type sensor. Then measured signal can be obtained by a photosensitive element, information acquired through signal will be shown on two groups of LEDs through process in electronic circuits and microprocessor.

Diagram of Operation Principle

1	2	3

1. One hang lace 2. One user manual

Guidance and manufacturer's declaration –electromagnetic emission –for all EQUIPMENT AND SYSTEMS

1	Guidance and manufacturer's declaration –electromagnetic immunity
2	The Fingertip Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the Fingertip Pulse Oximeter should assure that it is used in such an environment.
3	Immunity tests
4	RF emissions CISPR 11
5	RF emissions CISPR 11
6	Harmonic emissions IEC 61000-3-2
7	Voltage Fluctuations/Ricker emissions IEC 61000-3-3


Guidance and manufacturer's declaration –electromagnetic immunity –for all EQUIPMENT AND SYSTEMS

Guidance and manufacturer's declaration –electromagnetic immunity	
The Fingertip Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the Fingertip Pulse Oximeter should assure that it is used in such an environment.	
Immunity test	EN 60601 test level
Electrostatic discharge IEC 61000-4-3	2 kV contact
Electrostatic transient / burst IEC 61000-4-4	2 kV for power supply lines 2 kV for I/O lines
Surge IEC 61000-4-5	1 kV differential mode 2 kV common mode
Voltage dips, voltage variations on short interruption of power supply IEC 61000-4-11	5 kV for 5 cycles 5 kV for 5 cycles 5 kV for 5 cycles 5 kV for 5 cycles
Power frequency magnetic field IEC 61000-4-8	30 A/m

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
Precautions for use

1. Do not use the Fingertip Pulse Oximeter together with MRI or CT equipment.
2. Explosion hazard: Do not use the Fingertip Pulse Oximeter in an explosive atmosphere.
3. The Fingertip Pulse Oximeter is intended only as an adjunct in patient assessment. Doctors should make diagnosis in conjunction with clinical manifestation and symptoms.
4. The Fingertip Pulse Oximeter sensor application site frequently to make sure that the circulation and skin integrity of patient are under good condition.
5. Do not match the adhesive tape while applying the Fingertip Pulse Oximeter sensor. This may cause inaccurate reading or skin blisters.
6. Please read the manual carefully before your operation.
7. The Fingertip Pulse Oximeter has no SpO2 alarm, it is not for continuous monitoring.
8. When used for a long time, it will produce discomfort or tenderness, especially for patients with microcirculation disorders. Please change the test sites periodically according to different situations of patients. The test site must be changed and the skin integrity and circulation condition of the patient must be checked at least every 2 hours, and the correct adjustment must be made.
9. Inaccurate measurement may be achieved by autoclaving, ethylene oxide sterilizing, or immersing the sensors in liquid.
10. Significant levels of dysfunctional hemoglobins (such as carboxyl- or methemoglobin) may cause inaccurate reading.
11. Intravascular dyes such as indocyanine green or methylene blue may cause inaccurate reading.
12. SpO2 measurements may be adversely affected in the presence of high ambient light. Please shield the sensor area (with a surgical towel or direct sunlight, for example) if it is necessary.
13. Unexpected action may cause inaccurate reading if inserted into the finger, instrument is automatically switched on timely).
14. Medical signal with high frequency or interference caused by defibrillator may cause inaccurate reading.
15. Venous pulsations may cause inaccurate reading.
16. It may cause inaccurate reading when the positions of sensor and blood pressure cuff are not correct.
17. Hypotension, severe vasoconstriction, severe anemia, or hypothermia may cause inaccurate reading.
18. It may cause inaccurate reading by giving use of cardiotoxic to patient after his cardiac arrest or when he is in quiver.
19. Bright nail or painted nail may cause inaccurate SpO2 reading.

Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.

2

Guidance and manufacturer's declaration –electromagnetic immunity –for EQUIPMENT and SYSTEM that are not LIFE-SUPPORTING

Guidance and manufacturer's declaration –electromagnetic immunity	
The Fingertip Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the Fingertip Pulse Oximeter should assure that it is used in such an environment.	
Immunity test	EN 60601 test level
Conducted RF IEC 61000-4-6	150 VHz to 80 kHz
Compliance level	N/A

Table 9 – Test specifications for enclosure port immunity to RF wireless communications equipment

Test frequency (MHz)	Bandwidth (MHz)	Test wave form	Modulation (b)	Maximum power (W)	Distance (m)	Immunity level
385 – 390 – 390	TETRA 400	Pulse modulation (b)	1.8	0.3	27	2
450 – 430 – 470	GSM 400, FRS 400	FM/CT ± 5 kHz deviation	2	0.3	28	2
746 – 704 – 787	LTE Band 13,17	Pulse modulation (b)	2	0.3	28	2
800 – 800 – 900	GSM 800/900,TETRA 800, DECT 800,CDMA 800, 1.25 MHz	Pulse modulation (b)	2	0.3	28	2
1710 – 1710 – 1900	GSM 1800,CDMA 1800, GSM 1900,CDMA 1900, 1.25 MHz	Pulse modulation (b)	2	0.3	28	2
2450 – 2400-510	Bluetooth,WLAN 802.11 b/g/n, RFID 2450, LTE Band 1	Pulse modulation (b)	2	0.3	28	2
2400 – 2400 – 500	WLAN 802.11 a/n	Pulse modulation (b)	0.2	0.3	9	2

NOTE: If necessary to achieve the immunity test level, the distance between the transmitting antenna and the me equipment or me system may be reduced to 0.1 m. The 1 m distance is permitted by IEC 61000-4-3. a) For some systems, only the uplink frequencies are included. b) For some systems, the test may be conducted with the modulated signal. c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

Possible Problems and Resolutions

Problem	Possible reason	Solution
SpO2 or PR can not be shown normally	1. Finger is not plugged correctly 2. Patient's hemoglobin value is too low to be measured	1. Retry by plugging the finger 2. Try more times. If you can make sure there is no problem in product, please go to hospital timely for exact diagnosis
SpO2 or PR is shown unstable	1. The Finger might not be plugged deep enough 2. Finger is trembling or the patient is on movement status	1. Retry by plugging the finger 2. Please remain at rest

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Scope of application / Intended use

The fingertip pulse oximeter can be used to measure human haemoglobin saturation and pulse rate through finger , it can be used in hospitals, families, schools and medical centers.

contraindication:

note

- 1.The image in the instruction may have slight differences with the actual instruments.
- 2.Technical parameters and appearance change, without prior notice.

Product include: main machine and SPO2 sensor.

Features

- ◆OLED display
- ◆Low voltage indicator
- ◆If the hand movements, under the effect of accelerometer, the interface can have four different kinds of display mode (suitable for medical accelerometer function instrument).
- ◆Low-power consumption, continuously work for more than six hours with two AAA batteries
- ◆In the absence of signals, the product will be in after 8 seconds to enter standby state
- ◆Small in volume, light in weight, and convenient to carry
- ◆Instrument has 5s automatic signal detection function, when you insert finger, timely will automatically start;Automatic startup function instrument (applicable to Automatic startup function instrument)

Operation Instructions

- 1.Install two AAA batteries into battery cassette before covering its cover.
- 2.Plug one finger into rubber hole of the Oximeter (It is best to plug the finger thoroughly) before releasing the clamp with the nail upwards.
- 3.Press button on the front panel; (Note: if equipped with automatic startup function refers to clamp oximeter; need not press the button, the instrument has 5 s automatic signal detection function, when you insert the finger, instrument is automatically switched on timely).
- 4.SpO2 and PR, respectively, will be displayed on the screen.
- 4.Don't tremble your finger when the Oximeter is working. Your body is not recommended on moving status.
- 5.Press the button on the front panel, if we want change display direction.
- 6.If equipped with accelerometer function of instrument then don't press the button, hand movements, the instrument with the accelerometer has four corresponding interface switch)
- 6.Read relevant datum from the display screen.
7. The instrument has the function of sleep, no signal will enter standby state of sleep;
- 8.Please replace new batteries when OLED indicates the batteries are in low power.

The Oximeter can not be turned on	1. Inadequate power or power off	1. Please replace the batteries
	2. Batteries might be installed incorrectly	2. Replace the batteries
	3. The Oximeter might be damaged	1. Normal
Indication lamps are suddenly off	1. The product automatically shuts off when no signal is detected for 8 seconds	2. Replace the batteries
	2. Inadequate power	

Reserves the right to technical change appearance, our products are subject to change without prior notice, please forgive me!

Statement:

Consentment with data such as circuit diagram, components list, figure and the detailed rules of correction, injection, available only to the repair factory training qualified personnel and units. The company can be in the form of email or other electronic files provide users with random files.

3. The instrument is not used for evaluation of blood oxygen probe pulse and pulse blood oxygen monitor accuracy.

After-sales service

Ensure that users

- Please read user manual before using the instrument;
- According to the requirement of the instruction manual for the operation and daily maintenance, and make sure the machine power supply, and environmental requirements

Maintenance time

Monday to Friday 9 to 17:30PM, except for the national legal holiday

Maintenance regulations

- To conform to the regulations, free maintenance within the scope of products, with warranty card for free maintenance.All that is beyond the scope of free maintenance product, provide paid services.
- With warranty card and shopping invoice,main machine for a year, accessories for three months are under free maintenance services from the date of purchase.
- Following does not belong to the scope of free maintenance
« The fault caused by human factors, the damage;
« Due to the use to be inconsistent with the provisions of our company work environment to cause damage;
« Due to the product in the our company authorized personnel disassembling or repairing damages;
« Products beyond the warranty period.
If any problem, please call us in 9: 00 am to 5: 30 pm from Monday to Friday (national holiday), call us:
+86-516-68910112 +86

Settings	
Alarm Setup	+
Alarm	off
Bleep	off
Demo	off
Restore	OK
Brightness	OK
Exit	3

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