Bilancia pesapersona digitale

Manuale d'Uso

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 Zhongshan Jinli Electronic Weighing Equipment Co., Ltd
Zhongshan City, Guangdong, China

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Digital Weight Scale

User Manual

This weight scale allows you to evaluate your body weight. The scale capacity is 150 kilos with an accuracy of 100 grams.

-  **Please read through these instructions carefully before using the scale.**

RECOMMENDATIONS FOR USE

- Place the scale on a flat, hard surface. Soft, uneven flooring (e.g. rugs, linoleum) can cause the scale to read the weight incorrectly.
- Step onto the scale with your feet parallel and weight equally distributed.
 - Stand still while the scale performs the measurement.
 - Always weigh yourself on the same scale each day at the same time, preferably undressed and before breakfast. To get the best results from your scale, weigh yourself twice, and if the two weights are different from each other, calculate the average of the two readings.
 - The surface of the scale may be slippery, if wet.

USING THE SCALE FOR THE FIRST TIME

How to insert the battery

This scale operates with one lithium battery. When the battery needs to be replaced the display shows «**Lo**» (1). Insert a new battery in the battery compartment (5) with the positive pole facing upward. Remove the battery if the scale is not used for a prolonged period.

Switching the Unit of Measure Kg/Lb/St

Located at the bottom of the scale, this switch (6) allows choosing the unit of measure for

weighing: Kg, Lb or St.

Select your desired unit of measure when the scale is turned on and «**0.0**» is displayed (otherwise the unit of measure cannot be changed).

Using the Scale

- To get accurate measurements, place the scale on a flat and stable floor.
- Step onto the scale and wait without moving; after a few seconds, the display screen will show your weight (3). Do not bend or move while being weighed as this could alter the reading.
- Step off the scale, the display will continue to show your weight for one minute then it will turn off automatically.
- If the «**Err**» symbol (4) is displayed, it means that an error has occurred during the measurement. Perform the weighing procedure again from the start.

SAFETY, CARE, ACCURACY TEST AND DISPOSAL

- This device may only be used for the purposes described in this manual. The manufacturer shall not be liable for any damage caused by incorrect use.
- This device comprises sensitive components and must be handled with care. Observe the storage and operating conditions described in the «Technical Specifications» section.
- Protect the scale from:
 - water and moisture
 - extreme temperatures
 - impact and falls
 - contamination and dust
 - direct sunlight
 - heat and cold
- Do not use this device if you think it is

EN

- damaged or notice anything unusual.
- Never open this device.
- Remove the battery if the device is not used for a prolonged period.

Device care

Clean the device only with a soft, dry cloth.

Disposal

 Batteries and electronic devices must be disposed of in accordance with  local applicable regulations, not with domestic waste.

TECHNICAL SPECIFICATIONS

Capacity:	150 Kg / 330 Lb / 24 St
Accuracy:	0.1 Kg / 0.2 Lb / 1/4 St
Power:	One 3V lithium battery, CR 2032
Tolerance:	+/- 1% + 0.1 Kg / 0.2 Lb / 1/4 St

The manufacturer reserves the right to apply technical changes.

WARRANTY

This device is covered by a 2-year warranty from the date of purchase, against manufacturing defects, valid only upon submission of the purchase receipt. Opening or tampering with the device voids the warranty. The warranty does not cover accidental damage caused by improper use or failure to abide by the instructions for use.

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