# Easy 2

PRO

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# Easy 2

PRO

Smart nebulizer  
synchronized with breathing

Read the instructions carefully before using this device.

Dear Customer, This nebulizer has been designed to treat chronic illnesses that require maximum medicine inhalation. It has an innovative breath-synchronised valve system. It can be used with all compressed air aerosols and it nebulizes all aerosoltherapy medicines. If you have any questions or problems or if you need to order spare parts, please contact your dealer or the COLPHARMA customer care service. Alternatively, you can visit [www.colpharma.com](http://www.colpharma.com) where you will find lots of useful information about our products.

## PREPARATION AND USAGE OF THIS DEVICE

When using the nebulizer or the nebulizer kit for the first time, make sure to clean all parts as described in the section «Cleaning and Disinfection».

1. Open the nebulizer and pour in the amount of medicine as prescribed by the physician. Make sure not to pass the maximum level mark (PIC. 1/A).
2. Choose the most suitable inhalation accessory, as prescribed by the physician. The breath-synchronised valve system maximises the effect of the medicine, and works only with the mouthpiece (PIC.1/B).
3. When using the mask or to shorten the treatment time with the mouthpiece, Open the cap with the inlet valve. The nebulizer will work with continued flow and therefore, more quickly (PIC. 1/C).
4. Connect the nebulizer to the compressor using the air tube (PIC. 1/D).
5. Press the ON/OFF button to switch on the compressor. Place the mouthpiece (3) in the mouth,

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rest the mask (8) against the face, making sure that it covers nose and mouth.

6. Sit in a relaxed position with the upper body straight. Do not lie down while inhaling.
7. Make sure that the air tube contains no condensation. In case of condensation, switch on the compressor to dry up any moisture in the tube and then disconnect it as described above.
8. Empty out the remaining medicine from the nebulizer and clean all the parts, as described in the section «Cleaning and Disinfection».

## CLEANING AND DISINFECTION

Carefully clean all parts and remove any medicine residues and possible impurities after each treatment. The air tube (9) must be wiped clean with a clean damp cloth.

Always wash hands thoroughly before cleaning and disinfecting the accessories.

Wash/sterilize every part of the nebulizer.

Replace the air tube when passing to a new patient or in case of impurities.

### Cleaning with water

- Clean all the parts (except for the air tube) under hot running water for about 5 minutes, adding a small amount of mild detergent if necessary, following the dosage and limits specified by the detergent manufacturer.
- Rinse thoroughly, taking care to remove all residues, then leave to dry.

### Disinfection

- All parts (except the air tube) can be disinfected with chemical disinfectants, using doses and limits as specified by the manufacturer of the disinfectant.
- Disinfectants are generally available from pharmacies.

### Steam sterilization

- All parts (except the air tube and the face masks) can be sterilized using steam of up to 100 °C (15-20 min.).
- After sterilizing, always leave the components to cool down to ambient temperature before using again.

Do not repeat the sterilisation cycle when the components are still hot.

### Sterilization in boiling water

- All of the components (except for the air tube and the

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masks) can be immersed in boiling water for no longer than 1 minute.

Boil the components in plenty of water.

Disassemble all components of the nebulizer and boil them separately.

## MALFUNCTIONS AND WHAT TO DO

- Make sure that the compressor is working.
- Make sure that the air tube (9) is correctly connected at both ends.
- Make sure the air tube is not squashed, bent, soiled or blocked. If necessary, replace with a new one.
- Make sure the nebulizer is completely assembled (PIC. 1/B/C) and ensure the components are correctly assembled.
- Make sure the medicine has been poured into the nebulizer.

## TECHNICAL DATA

Max filling volume: 8 ml  
Min filling volume: 2 ml  
Residual volume: 0,5 ml

## GUARANTEE

- The spare parts such as nebulizer, masks, mouthpiece, and air tube are guaranteed for one year for defects in materials and workmanship, but not against deterioration caused by normal wear and tear.
- The guarantee does not cover damage caused by improper use, accidents, or failure to abide by the instructions for use.

\* **The masks for adults and kids and the air tube are supplied only with the EASY 2 PRO KIT (code 10758) and not with the EASY 2 PRO (code 10760) nebulizer.**

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- EN**
- 1 Cap and intake valve
  - 2 Nebulizing system
  - 3 Mouthpiece with expiratory valve
  - 4 Accessories fitting
  - 5 Nebulizer body
  - 6 Nebulizer support / carry handle
  - 7 Air tube fitting
  - 8 Adult and Child masks\*
  - 9 Air tube\*

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